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RNA interference in mammalian cells using low siRNA concentrations

Jörg Dennig*, Silvia Magyar*, Anja Grewe*, Cornelia Schmidt*, Peter Hahn*, Dong Liang†, Subu Yerramilli†, Eric Lader†, Wolfgang Bielke*, and Jie Kang*.

* QIAGEN GmbH, Hilden, Germany † QIAGEN Sciences, Germantown, MD, USA


Introduction

The use of short interfering RNA (siRNA) for knockdown of gene expression has become a powerful tool in molecular and cell biology. Some applications require the use of low siRNA concentrations (less than 5 nM); for example, to decrease the possibility of nonspecific effects.

We have developed a transfection reagent, HiPerfect Transfection Reagent, which allows efficient gene knockdown with siRNA concentrations from 1 nM–10 nM, depending on the cell type and siRNA used. HiPerfect Transfection Reagent has been tested and validated for many cell types, including primary cells. Effective knockdown in primary cells demonstrates that HiPerfect Transfection Reagent ensures low cytotoxicity levels. A FastForward siRNA Transfection Protocol has been developed for rapid transfection with HiPerfect Transfection Reagent. This protocol allows cell seeding and transfection on the same day.

A reverse transfection protocol has been developed that is ideal for use in high-throughput applications. In reverse transfection, siRNA is spotted into wells, followed by addition of HiPerfect Reagent. After complex formation, cells are added to the wells.

Reverse Transfection Using HiPerfect Transfection Reagent

Highly effective knockdown of CDC2 expression with low siRNA concentrations

Comparison of Knockdown Efficiency Using HiPerfect Transfection Reagent and Reagent L


Figure 1 Hela 50 cells were transfected with a range of concentrations of siRNA targeted against CDC2 using HiPerfect Transfection Reagent from QIAGEN or Reagent L from another supplier. Knockdown efficiency (measured by qPCR) was also tested. After 48 hours, gene silencing was measured by qPCR. siRNA concentrations of 1 nM–10 nM were tested. For CDC2 analysis, total cellular RNA was prepared using the RNeasy Lysis Buffer, and total RNA was extracted using RNeasy Spin Columns. The resulting RNA was used for real-time RT-PCR. Knockdown of CDC2 expression was measured by qPCR. The results show that HiPerfect Transfection Reagent from QIAGEN allowed highly efficient CDC2 knockdown with siRNA concentrations as low as 1 nM.

HiPerfect Transfection Reagent


Figure 2 HiPerfect Transfection Reagent from QIAGEN allowed highly efficient CDC2 knockdown with siRNA concentrations as low as 1 nM. In contrast, Reagent L from another supplier provided less efficient knockdown of all concentrations tested. For concentrations lower than 5 nM, knockdown of only 15% or less was observed.

Transfection and knockdown in a wide range of cell types

A wide range of cell types have been successfully transfected using HiPerfect Transfection Reagent. For an up-to-date list of cell types and more detailed information go to www.qiagen.com/TransfectionTools.

A Wide Range of Successfully Transfected Cells


Figure 3 The complete cell list is available at www.qiagen.com/TransfectionTools.

Rapid, efficient lamin A/C knockdown in human primary cells

Lamin A/C Knockdown Using the FastForward Protocol with Low siRNA Concentrations


Figure 3 Human foreskin fibroblasts (HFFs) were transfected with siRNA targeted against human A/C using HiPerfect Transfection Reagent with the FastForward Protocol. The knockdown efficiency was analyzed by qPCR. siRNA concentrations of 1 nM–10 nM were tested. For A/C analysis, total cellular RNA was prepared using the RNeasy Lysis Buffer, and total RNA was extracted using RNeasy Spin Columns. The resulting RNA was used for real-time RT-PCR. Knockdown of A/C expression was measured by qPCR. The results show that HiPerfect Transfection Reagent from QIAGEN allowed highly efficient lamin A/C knockdown with siRNA concentrations as low as 1 nM.

Microscopic Analysis of Transfected Cells

Figure 4 Phase contrast microscopy of transfected cells.

- HiPerfect Transfection Reagent allowed highly efficient lamin A/C knockdown with siRNA concentrations as low as 1 nM.
- The FastForward Transfection Protocol allowed rapid transfection, with cell seeding and transfection carried out on the same day.
- Low cytotoxicity means that HiPerfect Transfection Reagent is especially suitable for use with sensitive primary cells.

HiPerfect Transfection Reagent Allows Effective Uptake of Low Amounts of Alexa Fluor® 488 Labeled siRNA

Fluorescently Labeled siRNA Shows Highly Efficient Uptake at Low siRNA Concentrations


Figure 5 Fluorescence microscopy of HeLa cells 24 hours after transfection with 2 nM Alexa Fluor 488 labeled siRNA (green fluorescent) using HiPerfect Transfection Reagent. Inside most stained cells nucleoli (blue) are visible.

Summary

- HiPerfect Transfection Reagent allows gene silencing using as little as 1 nM siRNA. Transfection of low siRNA concentrations may be necessary to avoid off-target effects. Using HiPerfect Transfection Reagent means that effective knockdown can be achieved with very low siRNA concentrations.
- HiPerfect Transfection Reagent has been tested for a range of cell types. Many cell lines and primary cells have been successfully transfected using low siRNA concentrations and HiPerfect Transfection Reagent. For an up-to-date list, visit www.qiagen.com/TransfectionTools.
- The FastForward Protocol has been developed for rapid transfection. In the FastForward Protocol, cell seeding and siRNA transfection are carried out on the same day, saving time and effort. The FastForward Protocol is available at www.qiagen.com/qiagen/HiPerfect.
- The reverse transfection protocol is ideal for RNAi screening. The reverse transfection protocol can be easily automated which is particularly useful for high-throughput applications. The reverse transfection protocol is available at www.qiagen.com/qiagen/HiPerfect.

Transfection: QIAGEN's QIAGEN HiPerfect Transfection Reagent allows HiPerfect Transfection Reagent. The HiPerfect protocol is covered by U.S. Patents 6,855,151 and 6,855,152 and foreign equivalents owned by QIAGEN. QIAGEN's HiPerfect technology based in Germany is covered by various patent applications, owned by the QIAGEN family of companies. QIAGEN, Hilden, Germany.

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

PROBITY

IAN W. DAVIS, W. BRYAN ARENDALL III, LAURA W. MURRAY, JEREMY N. BLOCK,
JANE S. RICHARDSON, DAVID C. RICHARDSON
DEPARTMENT OF BIOCHEMISTRY, DUKE UNIVERSITY, DURHAM, NC

THINGS THAT GO "BUMP" IN PROTEINS:

...AND HOW TO FIND THEM:


OVERVIEW


The PROBITY pipeline identifies and tracks the most common conformational changes in proteins, and is designed to be used by researchers in a variety of contexts. It can be used to identify conformational changes in a single protein, or to compare conformational changes across a set of proteins. It can also be used to identify conformational changes in a specific region of a protein, or to identify conformational changes in a specific type of protein. The pipeline is designed to be used by researchers in a variety of contexts, and is designed to be used by researchers in a variety of contexts.

NEW! FASTER HYDROGENS!

Residue Distributions for Residue Variations

RESULTS @ SEC66

RESULTS WITH RNA

[HTTP://KINEMAGE.BIOCHEM.DUKE.EDU](http://KINEMAGE.BIOCHEM.DUKE.EDU)

CONCLUSIONS & REFERENCES

CONCLUSIONS: The PROBITY pipeline identifies and tracks the most common conformational changes in proteins, and is designed to be used by researchers in a variety of contexts. It can be used to identify conformational changes in a single protein, or to compare conformational changes across a set of proteins. It can also be used to identify conformational changes in a specific region of a protein, or to identify conformational changes in a specific type of protein. The pipeline is designed to be used by researchers in a variety of contexts, and is designed to be used by researchers in a variety of contexts.

REFERENCES: Davis, I. W., Arendall, W. B., Murray, L. W., Block, J. N., Richardson, J. S., Richardson, D. C. (2004) PROBITY: A pipeline for identifying and tracking the most common conformational changes in proteins. *Journal of Molecular Biology*, 341, 1-10.

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